



Anant, A. Reka
Photography




Sample characterization using X-ray powder diffraction technique

Assoc. prof. dr. Arianit A. Reka

¹Department of Chemistry, Faculty of Natural Sciences and Mathematics, University of Tetova, Tetovo, 1200, North Macedonia

²NanoAlb - Unit of Albanian Nanoscience and Nanotechnology, Albanian Academy of Sciences, Tirana, Albania



23 march 2023

Content:

- Why XRPD?
- How does it work?
- Interpretation...
- Questions?
- Thank you note

- X-ray powder diffraction is most widely used for the identification of unknown crystalline materials (e.g. minerals, inorganic compounds) in environmental science, material science, engineering, biology, etc.




Figure 1. Clayey diatomite, pure diatomite and non-expanding perlite, fly-ash, bentonite clay, magnesite

How?








Figure 2. XRPD diffractometer

How?



Amorphous Solids



Crystalline Solids

Figure 3. Amorphous vs crystalline (constructive interference)

Strengths

- Powerful and rapid (< 20 min) technique for identification
- Minimal sample preparation is required
- XRD units are widely available
- Data interpretation is relatively straight forward
- Non-destructive technique

Limitations

- Sometimes it requires tenths of a grams of material which must be ground into a powder
- For mixed materials, detection limit is ~ 2-5% of sample
- Must have access to a standard reference file of inorganic compounds (d-spacings)





Figure 4. Sample of hematite ore



<https://handbookofmineralogy.org/>

<http://webmineral.com/>

<https://rruff.info/>

Portlandite




©2001-2005 Mineral Data Publishing, version 1

X-ray Powder Pattern: Synthetic.

2.628 (100), 4.90 (74), 1.927 (42), 1.796 (36), 3.112 (23), 1.687 (21), 1.484 (13)

Quartz



©2001 Mineral Data Publishing, version 1.2

X-ray Powder Pattern: Synthetic.


3.342 (100), 4.257 (22), 1.8179 (14), 1.5418 (9), 2.457 (8), 2.282 (8), 1.3718 (8)





SUSNANO

THANK YOU!



NanoAlb



Funded by the
European Union

